EN INFORMATION FOR USER

VIZOOVET® PROTECT™
Sterile Eye Drops 6 ml
10 Monodose vials 0,6 ml
Relief and freshness for irritated eyes

INDICATIONS: Eye drops for stressed, irritated or tired eyes, intended to wash and hydrate eyes as physiological response to external factors.
Composition: Purified water, glycerin, propolis extract aloe vera extract, matricaria chamomile extract, sodium chloride, boric acid, sodium tetraborate, EDTA disodium.

RECOMMENDATIONS FOR USE: Apply 1-2 drops to each eye. Opened vial is recommended for single use. Throw away in 24 hours after vial opening.
Safe to use as often as needed.
Do not use on bandage lenses directly; wait at least 10 minutes after using drops before inserting bandage lenses.

PRECAUTIONS:
Do not touch the dropper of vial on to any surface to avoid contamination.
Do not ingest and keep out of the reach of children.
Do not use if you are allergic to any of the components.
Do not store above 30° C.
Do not freeze.
Do not use if there are visibles particles in the solution.
Do not use if vial is damaged or leaking.
Do not use after the expiry date which is printed on the carton and vial after EXP.
The expiry date refers to the last day of the month.
In case of irritation consult your veterinarian.
Follow all directions on the product package and information for user or use as directed by your veterinarian.

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VIZOOVET® PROTECT™
Gocce oculari sterili da 6 ml
10 flaconi monouso da 0,6 ml
Rinfresca e da sollievo agli occhi irritati

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